

Assay Kit versions

All Assay Kits include 25 Biosensors, ready-to-use Buffer and Positive Control

	LactoSens®R	LactoSens®R for NOLA™ Fit	LactoSens® 0.02 %
Measurement range	0.008 - 0.2 % lactose	0.008 - 0.2 % lactose	0.02 - 1.0 % lactose
Lactase	Yeast neutral (e.g.: HA-Lactase, Opti-Lactase LX2, Maxilact)	Bifidobacteria (e.g.: NOLA™ Fit, Eco-B 5000, Saphera)	Yeast neutral (e.g. HA-Lactase, Opti-Lactase LX2, Maxilact)
Sample dilution	1+1	1+1	1+4
NordVal Certification	yes	in progress	no
Application notes*			
Milk/Cream	yes	yes	yes
Flavored milk	yes	yes	no
Yogurt/Sour cream	yes	yes	yes
Flavored/Fruit Yogurt	yes	yes	no
Cottage/Curd/Fresh cheese	yes	in progress	yes
Mozzarella/Hard cheese	yes	in progress	yes

* Register on website for download or contact your distributor