

Compare Assay Kit versions

All Assay Kits include 25 Biosensors, ready to use dilution buffer and positive control

	LactoSens® R	LactoSens® R for NOLA™ Fit	LactoSens® 0.02 %	LactoSens® 0.01 %**
Measurement range	0.008 - 0.2 % lactose	0.008 - 0.2 % lactose	0.02 - 1.0 % lactose	0.01 - 0.2 % Lactose
Lactase	Yeast neutrol (e.g.: HA-Lactase, Opti-Lactase LX2, Maxilact)	Bifidobacteria (e.g.: NOLA Fit, Eco-B 5000, Saphera)	Yeast neutral (e.g. HA-Lactase, Opti-Lactase LX2, Maxilact)	Yeast neutral (e.g.: HA-Lactase, Opti-Lactase LX2, Maxilact)
Sample dilution	1+1	1+1	1+4	1+1
NordVal Certification	yes	In progress	no	yes
Application notes*				
Milk and Cream	yes	yes	yes	yes
Flavored milk	yes	yes	no	no
Yogurt and Sour cream	yes	yes	yes	yes
Flavored/Fruit Yogurt	yes	yes	no	no
Cottage, Curd and Fresh cheese	yes	In progress	yes	no
Mozzarella and Hard cheese	yes	In progress	yes	no

* Register on website for download or contact your distributor

** Product discontinued - available only until July 2019